



PINK Color Kinesiology Tape Medical Muscle Tape 5cmx5m Reduce Pain CE ISO Approved

Our Product Introduction

for more products please visit us on medicalcare-products.com

Basic Information

- Place of Origin: China
- Brand Name: JaneJoy
- Certification: FDA and CE
- Model Number: J-04
- Minimum Order Quantity: 1000rolls
- Price: 0.85
- Packaging Details: 6rolls in innerbox, 144rolls in one carton
- Delivery Time: 30days lead time
- Payment Terms: T/T
- Supply Ability: 100000000rolls per month




Product Specification

- Material: 95% Cotton+5%Spandex
- Size: 5cm*5m
- Color: Pink And 22colors For Choice
- Certificate: FDA And CE
- Type: Muscle Tape
- MOQ: 1000rolls
- Custom: Customization Offered
- Packing: 144rolls One Carton
- Highlight: **5cmx5m Medical Muscle Tape, ISO Approved Medical Kinesiology Tape, CE Approved Kinesiology Tape**



More Images



Product Description

PINK Color Kinesiology Tape Muscle Adhesive Medical Tape CE ISO Approved


Key attributes

Other attributes

Place of Origin	Jiangsu, China
Brand Name	BYTECH or FITOP
Model Number	5CM*5M
Function	Reduce Pain
Name	kinesiology tape
Material	95% Cotton+5%Spandex
Certificate	CE/ ISO
Type	Muscle tape
MOQ	144 rolls
Color	Pink
Delivery Time	7-15 Days
Sample	Free sample available
Custom	Customization Offered

Products Description





Product Information	
Product name	Kinesiology Tape/Athletic Tape/Sports tape
Material	95% Cotton+5%Spandex
Elasticity	Stretch to 170%-180%
Fabric weight	165 g/sqm

MOQ	120 rolls or 144 rolls, just one carton to be discussed
Peeling Strength	1.4 n/cm - 1.8 n/cm
Certification	CE/ISO9001/ISO13485
Adhesive Strength	5N
Size	2.5cm*5m,3.8cm*5m,5cm*5m,7.5cm*5m,10cm*5m or Customized
Color	12 colours are available or Customized
Logo	Either can be available in your logo or our own logo,as your request.OEM
Free sample	Acceptable
Used for	Body , Arm , leg , foot , neck , back , waist etc.
Features	High Sticky,waterproof, breathable and durable
Glue	Skin Friendly Acrylic Medical Glue
Function	Promote blood circulation,reduce muscle tension,reduce pain
Usage	Training, Fitness, Sports,Medical beauty
Quality	Each product will be 100%checked before shipment
Trade Information	
Sample order	About 3-7days for prepare, and shipped by Express.
Bulk order	About 7-15 days , that's depends on the inventory, quantity and custom logo
Payment terms	Accept T/T,L/C, Western Union, MoneyGram,Paypal, ,Alibaba Trade Assurance
Shipping terms	By Express, air way, or sea freight, according to customer's requirement.

Packing & Delivery

--



Guangzhou Janejoy Medical Technology Co.,Ltd



+8613535187404



JaneJoy@therapy.org.cn



medicalcare-products.com

Building A2, Loteam Creative Park, No. 251, Kehua Street, Tianhe District, Guangzhou, China