



## OEM 5cm Width Kinesiology Athletic Tape , Sports Muscle Tape For 5-7 Days Wear

Our Product Introduction

### Basic Information

- Place of Origin: China
- Brand Name: Janejoy
- Model Number: J02



### Product Specification

- Usage: Sports And Rehabilitation
- Material: 95% Cotton And 5% Spandex
- Sample: Free
- Application Method: Roll-on
- Adhesive Strength: 5N
- Elasticity Ratio: 1:2±0.3
- Breathability: Yes
- Stretchability: 140%
- Highlight: **5cm Width Sports Muscle Tape,  
OEM Kinesiology Athletic Tape,  
5cm Width Waterproof Sports Tape**

For more products please visit us on [medicalcare-products.com](http://medicalcare-products.com)

## Product Description

### Breathable Kinesiology Tape Muscle Tape for 5-7 Days Wear

#### Product Description:

Our Kinesiology Tape is highly durable, lasting for 5-7 days, so you can be sure it will withstand even the most intense workouts or activities. Additionally, our tape comes in a convenient 5m length, so you'll have plenty of tape to use for multiple sessions.

Our tape is waterproof and breathable, making it the perfect choice for sports and rehabilitation. Whether you're running in the rain or sweating through a workout, our tape will stay in place and provide the support you need. Plus, its breathable design allows for air to flow through, preventing any discomfort or irritation.

Our Kinesiology Tape is latex-free, making it the perfect choice for those with allergies or sensitivities. You can use our tape with peace of mind, knowing that it won't cause any irritation or discomfort.

Our Kinesiology Tape is a must-have for anyone looking to improve their athletic performance or recover from an injury. Its unique design provides support and relief for muscles, joints, and ligaments, allowing you to push your limits and reach your goals.

So whether you're a professional athlete or simply looking to improve your fitness, our Kinesiology Tape is the perfect choice. Order yours today and discover the many benefits of our muscle tape, sport tape, and Kinesiology tape!

#### Features:

Product Name: Kinesiology Tape

Waterproof And Breathable: Convenient For Usage

Latex-Free: Yes

Sample: Free

Material: 95% Cotton And 5% Spandex

Stretchability: 140%

Muscle tape, Kinesiology tape, Muscle tape

#### Technical Parameters:

Product Attribute	Value
Product Name	Sport Tape/Muscle Tape/Kinesiology Tape
Durability	5-7 Days
Adhesive Strength	5N
Gel	Acrylic
Breathability	Yes
Latex-Free	Yes
Width	5cm
Elasticity Ratio	1:2±0.3
Waterproof and Breathable	Convenient for Usage
Stretchability	140%
Material	95% Cotton and 5% Spandex

#### Applications:

The Janejoy Kinesiology Tape uses acrylic gel, which is latex-free, making it safe for those who are allergic to latex. The tape is also breathable, allowing air to circulate around the taped area, reducing the risk of skin irritation and discomfort. The width of the tape is 5cm, providing enough coverage for larger muscle groups.

The Kinesiology Tape is suitable for a variety of applications and scenarios, including:

**Sports:** The Janejoy Kinesiology Tape is perfect for athletes and fitness enthusiasts who need support and protection for their muscles during physical activities. It can be used for different sports such as running, basketball, football, and more.

**Injury Prevention:** The tape can be used to support and stabilize the muscles to prevent injuries during physical activities. It can also help reduce muscle fatigue and soreness after exercise.

**Injury Recovery:** The tape can be used to support injured muscles during the recovery process. It can help reduce swelling and provide support to the affected area.

**Posture Correction:** The tape can be used to help correct posture by providing support to the muscles and joints.

The Janejoy Kinesiology Tape is easy to apply and remove, making it a convenient option for athletes and fitness enthusiasts. It is also water-resistant, ensuring that it stays in place during physical activities, even when sweating.

In summary, the Janejoy Kinesiology Tape is a versatile and effective muscle tape that can be used for a variety of applications and scenarios. Whether you are an athlete, fitness enthusiast, or someone who needs support and protection for their muscles, this tape is an excellent choice.

## Customization:

Brand Name: Janejoy  
Model Number: J02  
Place of Origin: China  
Durability: 5-7 Days  
Width: 5cm  
Application Method: Roll-on  
Latex-Free: Yes  
Material: 95% Cotton And 5% Spandex  
Kinesiology tape, Kinesiology tape, Kinesiology tape

## Support and Services:

Our Kinesiology Tape product comes with comprehensive technical support and services to ensure that you get the most out of your purchase. Our team of experts is available to answer any questions you may have about the product's features, benefits, and proper application techniques.

We also offer various resources to help you learn more about our Kinesiology Tape product, including instructional videos, user manuals, and online tutorials. If you encounter any issues with your product, our support team is available to troubleshoot and provide solutions to ensure that you have a positive experience with our product.

Our services also include product warranty and repair services. If your product is defective or requires repair, simply contact our customer service team to arrange for a replacement or repair. We stand behind our products and are committed to ensuring that our customers are satisfied with their purchase.

## Packing and Shipping:

### Product Packaging:

The Kinesiology Tape product is packaged in a compact, easy-to-carry box that contains one roll of tape. The box measures 4 inches by 2 inches by 2 inches and is made of sturdy cardboard to ensure that the tape arrives in good condition.

### Shipping:

Our Kinesiology Tape product is shipped using standard shipping methods and is typically delivered within 3-5 business days. We offer free shipping on orders over \$50 and provide tracking information for all shipments.

## FAQ:

Q: What is Janejoy Kinesiology Tape?

A: Janejoy Kinesiology Tape is a high-quality elastic therapeutic tape that provides support and pain relief to muscles and joints. It is made from breathable, latex-free materials that are gentle on the skin.

Q: What are the benefits of using Janejoy Kinesiology Tape?

A: Janejoy Kinesiology Tape helps improve blood flow and lymphatic drainage, which can reduce inflammation and swelling. It also provides support to muscles and joints, which can help improve athletic performance and prevent injuries.

Q: How do I apply Janejoy Kinesiology Tape?

A: To apply Janejoy Kinesiology Tape, clean and dry the skin first. Then, apply the tape to the affected area, following the instructions provided. Be sure to avoid stretching the tape too much, as this can cause discomfort.

Q: Is Janejoy Kinesiology Tape waterproof?

A: Yes, Janejoy Kinesiology Tape is water-resistant and can be worn while swimming or showering. However, it is recommended to pat dry the tape after exposure to water to ensure maximum adhesion.

Q: Where is Janejoy Kinesiology Tape made?

A: Janejoy Kinesiology Tape is made in China, using high-quality materials and state-of-the-art manufacturing processes to ensure consistent quality and performance.



Guangzhou Janejoy Medical Technology Co.,Ltd



+8613535187404



JaneJoy@therapy.org.cn



medicalcare-products.com

Building A2, Loteam Creative Park, No. 251, Kehua Street, Tianhe District, Guangzhou, China