



## Breathable Kinesiology Pain Relief Muscle Tape Bandage 5cmx5m For Enhanced And Recovery

Our Product Introduction

### Basic Information

- Place of Origin: China
- Brand Name: Janejoy
- Model Number: J02



### Product Specification

- Elasticity Ratio: 1:2±0.3
- Adhesive Strength: 5N
- Length: 5m
- Width: 5cm
- Application Method: Roll-on
- Breathability: Yes
- Stretchability: 140%
- Latex-Free: Yes
- Highlight: 5cmx5m Pain Relief Muscle Tape Bandage,  
5N Adhesive Strength Muscle Tape,  
140% Stretchability Pain Relief Muscle Tape

For more products please visit us on [medicalcare-products.com](http://medicalcare-products.com)

## Product Description

### Breathable Pain Relief Muscle Tape Bandage For Enhanced Performance And Recovery

#### Product Description:

The Kinesiology Tape comes in a standard width of 5cm and a length of 5m, providing enough tape to cover large areas of the body. The width and length of the tape make it easy to apply and use on different parts of the body. This tape is perfect for individuals who need extra support and stability during physical activities, especially those involving the joints and muscles.

This muscle tape is made with high-quality materials that are latex-free, making it suitable for individuals with sensitive skin. The tape's adhesive is designed to be gentle on the skin, preventing irritation and discomfort. It is perfect for individuals who need a reliable and comfortable tape to use during physical activities.

If you are unsure whether the Kinesiology Tape is suitable for you, you can request a free sample to try it out. The free sample allows you to test the tape's strength, durability, and comfort before purchasing a full roll.

In conclusion, the Kinesiology Tape is a reliable and effective muscle tape that provides support and relief during physical activities. It is durable, easy to apply, latex-free, and comes with a free sample for individuals who want to test it out before purchasing. Whether you are an athlete or an individual who engages in physical activities, the Kinesiology Tape is an excellent choice for providing extra support and stability to your muscles and joints.

#### Features:

Product Name: Kinesiology Tape

Stretchability: 140%

Application Method: Roll-on

Latex-Free: Yes

Material: 95% Cotton And 5% Spandex

Length: 5m

This muscle tape, also known as kinesiology tape or sport tape, is designed with 140% stretchability and a roll-on application method. It is also latex-free and made of 95% cotton and 5% spandex, providing a comfortable and flexible fit. With a length of 5m, this kinesiology tape is perfect for athletes and active individuals seeking support for their muscles during physical activity.

#### Technical Parameters:

|                           |                           |
|---------------------------|---------------------------|
| Material                  | 95% Cotton And 5% Spandex |
| Stretchability            | 140%                      |
| Breathability             | Yes                       |
| Waterproof And Breathable | Convenient For Usage      |
| Width                     | 5cm                       |
| Length                    | 5m                        |
| Durability                | 5-7 Days                  |
| Application Method        | Roll-on                   |
| Latex-Free                | Yes                       |
| Gel                       | Acrylic                   |

#### Applications:

The adhesive strength of this kinesiology tape is 5N, which ensures that it stays in place during intense physical activities. Its elasticity ratio of 1:2±0.3 makes it ideal for providing support and compression to muscles and joints. The width of the tape is 5cm, which is wide enough to cover most areas of the body. Additionally, this sport tape is latex-free, making it safe for people with latex allergies to use.

This kinesiology tape is versatile and can be used in a variety of situations. It is commonly used to provide support and stability to muscles and joints during sports activities such as running, weightlifting, and basketball. It can also be used to alleviate pain and discomfort associated with injuries and illnesses such as arthritis, tennis elbow, and plantar fasciitis.

The Janejoy Kinesiology Tape is also suitable for post-surgery recovery. It can help reduce swelling, improve circulation, and promote faster healing. Additionally, it can be used for injury prevention by providing support to weak or vulnerable areas of the body.

In summary, Janejoy Kinesiology Tape is a high-quality muscle tape that offers support and stability to muscles and joints during physical activities. Its adhesive strength, elasticity ratio, and width make it suitable for a variety of situations. It is commonly used by athletes and fitness enthusiasts to prevent injury and alleviate pain. It is also suitable for post-surgery recovery and injury prevention. Made in China, this kinesiology tape is a reliable and effective tool for anyone looking to improve their physical performance and overall health.

#### Customization:

## Support and Services:

Our Kinesiology Tape product comes with comprehensive technical support and services to ensure that you get the most out of your purchase. Our team of experts is available to assist you with any questions or concerns you may have about the product, including:

Proper application techniques

Usage instructions

Product specifications

Troubleshooting common issues

In addition to our technical support, we also offer a range of services to help you get the most out of your Kinesiology Tape. These include:

Product training and education

Customized product solutions

Consulting and advisory services

Product customization and branding

We are committed to providing our customers with the highest level of support and service, and we stand behind our products with a satisfaction guarantee. Contact us today to learn more about our Kinesiology Tape technical support and services.

## Packing and Shipping:

Product Packaging:

The Kinesiology Tape comes in a durable plastic container with a secure lid to protect the tape from dust and moisture. The container is designed to be compact and easy to carry, making it convenient for athletes and individuals who need to use the tape while on the go.

Shipping:

We offer free shipping for orders of Kinesiology Tape within the United States. For international orders, shipping fees will be calculated based on the destination country and the weight of the package. All orders will be shipped within 24 hours of purchase and will be delivered within the estimated timeframe depending on the shipping method selected at checkout.

## FAQ:

**Q: What is Janejoy Kinesiology Tape?**

A: Janejoy Kinesiology Tape is a therapeutic tape made of elastic cotton cloth with an acrylic adhesive backing. It is designed to aid in the treatment of various physical injuries and conditions.

**Q: How does Janejoy Kinesiology Tape work?**

A: The tape works by lifting the skin and creating space between the skin and muscles. This helps to increase blood and lymph circulation, reduce inflammation, and alleviate pain and discomfort.

**Q: What conditions can Janejoy Kinesiology Tape be used for?**

A: The tape can be used to treat a variety of conditions, including muscle strains and sprains, joint pain, tendonitis, and swelling. It can also be used to improve athletic performance and prevent injuries.

**Q: Is Janejoy Kinesiology Tape easy to apply?**

A: Yes, the tape is easy to apply. Simply cut the desired length of tape, remove the backing, and apply it to the affected area. It can be worn for up to 5 days and is water-resistant.

**Q: Where is Janejoy Kinesiology Tape manufactured?**

A: Janejoy Kinesiology Tape is made in China.



Guangzhou Janejoy Medical Technology Co.,Ltd



+8613535187404



JaneJoy@therapy.org.cn



medicalcare-products.com

Building A2, Loteam Creative Park, No. 251, Kehua Street, Tianhe District, Guangzhou, China