



## Home Medical Portable Mesh Nebulizer Handheld Micro Grid Atomizer Rechargeable

Our Product Introduction

### Basic Information

- Place of Origin: China
- Brand Name: Janejoy
- Model Number: J04



### Product Specification

- Type: Ultrasonic
- Power Supply Mode: USB Recharge
- Customer Logo: Customized
- Warranty: 1 Year
- Instrument Classification: Class II
- Material: ABS
- Atomization Rate:  $\geq 0.2\text{ML/min}$
- Noise Level: Less Than 60dB
- Highlight: **Class II Handheld Micro Grid Atomizer, Customized Medical Portable Mesh Nebulizer, OEM Rechargeable Portable Mesh Nebulizer**

For more products please visit us on [medicalcare-products.com](http://medicalcare-products.com)

## Product Description

Nebulizer Children'S Home Medical Handheld Atomizer Adult Portable Micro-Grid Atomizer Rechargeable Atomizer

### Products Description

<b>Name</b>	Nebulizer
<b>Model</b>	U210
<b>Dimensions</b>	67*42*116mm
<b>Power</b>	0.75w
<b>Power supply</b>	USB charging
<b>Atomization rate</b>	>0.25ml/minute
<b>Atomized median particle size</b>	MMAD 4.0um
<b>Weight</b>	948g
<b>Ambient temperature</b>	10-40°C
<b>Maximum water temperature in medicine cup</b>	50°C Maximum
<b>Work noise</b>	50dB Maximum



Guangzhou Janejoy Medical Technology Co.,Ltd



+8613535187404



JaneJoy@therapy.org.cn



medicalcare-products.com

Building A2, Loteam Creative Park, No. 251, Kehua Street, Tianhe District, Guangzhou, China