



OEM Logo USB Powered Portable Mesh Nebulizer Machine For Respiratory Wellness

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: Janejoy
- Model Number: J04



Product Specification

- Material: ABS
- Type: Ultrasonic
- Instrument Classification: Class II
- Product Type: Portable Type
- Atomization Rate: $\geq 0.2\text{ML/min}$
- Power Supply: USB Cable
- Oem: YES
- Operate Method: Adapter For Power
- Highlight: **USB Powered Portable Mesh Nebulizer, OEM Logo Portable Nebulizer Machine, Portable Type Mesh Nebulizer Machine**

For more products please visit us on medicalcare-products.com

Product Description

OEM Convenient USB Powered Portable Mesh Nebulizer for Respiratory Wellness

Product Description:

Our Portable Nebulizer is powered by a USB cable, giving you the freedom to use it anywhere you go. The ultrasonic technology used in this nebulizer ensures that it is efficient and produces fine mist particles that can easily be inhaled.

One of the unique features of our Portable Nebulizer is that you can customize it with your own logo. This makes it an excellent option for businesses or healthcare providers who want to promote their brand while providing a useful product to their customers or patients.

The Portable Nebulizer is easy to use, with an adapter for power that allows you to plug it into any outlet. It is also incredibly lightweight and compact, making it easy to carry in a purse, backpack, or luggage.

We are proud to offer OEM services for our Portable Nebulizer, so you can customize it to fit your specific needs. Whether you need a different color, size, or design, we can work with you to create a nebulizer that meets your requirements.

Our Portable Nebulizer is also known as a Mesh Nebulizer, which is a type of nebulizer that uses mesh technology to create an aerosol mist. This type of nebulizer is preferred by many people because it is quiet, efficient, and produces a finer mist than traditional nebulizers.

If you're looking for a high-quality Portable Nebulizer that is easy to use and convenient to carry, then look no further than our product. Contact us today to learn more about our customization options and how we can help you with your respiratory therapy needs.

Features:

Product Name: Portable Nebulizer Machine

Atomization Rate: $\geq 0.2\text{ML/min}$

Product Type: Portable Type

OEM: YES

Type: Ultrasonic

Nebulizing Rate: $\geq 0.25\text{ml/min}$

Introducing the **Portable Nebulizer** - a convenient and efficient way to inhale medication. With an atomization rate of $\geq 0.2\text{ML/min}$ and nebulizing rate of $\geq 0.25\text{ml/min}$, this **portable nebulizer** is perfect for individuals on-the-go. Its ultrasonic technology ensures a fast and quiet treatment. Plus, it's OEM compatible for your convenience. Get the relief you need with the **Portable Nebulizer Machine**.

Technical Parameters:

Product Type	Portable Nebulizer
Instrument Classification	Class II
Material	ABS
Nebulizing Rate	$\geq 0.25\text{ml/min}$ (Mesh Nebulizer)
Noise Level	Less Than 60dB
Power Supply	USB Cable
Operate Method	Adapter For Power
Warranty	1 Year
OEM	YES
Customer Logo	Customized

Applications:

One of the standout features of the J04 Portable Nebulizer Machine is its use of mesh nebulizer technology. Unlike traditional nebulizers that use compressed air to break up liquid medication into a mist, mesh nebulizers use a vibrating mesh to create a fine aerosol. This results in a more efficient and consistent delivery of medication, with less waste and noise. The J04 Portable Nebulizer Machine's mesh nebulizer is also easy to clean and maintain, making it a great choice for busy individuals and families.

The J04 Portable Nebulizer Machine is perfect for a variety of application occasions and scenarios. It is an ideal choice for frequent travelers who need to maintain their respiratory health while away from home. Its compact size and lightweight design make it easy to pack in a suitcase or carry-on bag. It is also a great option for outdoor enthusiasts who want to enjoy activities such as hiking or camping without worrying about their respiratory condition.

The J04 Portable Nebulizer Machine is also a popular choice for busy professionals who need to manage their respiratory health while at work. Its quiet operation and discreet design make it easy to use in the office, without disturbing coworkers. It is also a great option for parents of children with respiratory conditions, as it allows them to administer medication quickly and easily, without causing discomfort or anxiety for their child.

The J04 Portable Nebulizer Machine is available for OEM customization, with the option to add a customer logo to the product. It operates using an adapter for power, making it easy to use in any location with an electrical outlet. Whether you're at home, at work, or on-the-go, the J04 Portable Nebulizer Machine is a reliable and convenient choice for respiratory therapy.

Customization:

Customize your Janejoy J04 Portable Ultrasonic Nebulizer Machine to your liking!

Product Attributes:

Brand Name: Janejoy

Model Number: J04

Place of Origin: China

Type: Ultrasonic

Customer Logo: Customized

Atomization Rate: $\geq 0.2\text{ML/min}$

Operate Method: Adapter For Power

Power Supply Mode: USB Recharge

With our customization services, you can now personalize your nebulizer with your own logo or design. Whether you're looking for a Mesh Nebulizer or a Portable Nebulizer, we've got you covered.

Support and Services:

The Portable Nebulizer Machine is a medical device designed to deliver medication directly to the lungs. It is easy to use and maintain, with a compact and lightweight design that makes it ideal for travel or home use. Our product technical support and services include:

Expert technical assistance for troubleshooting and problem-solving

Repair and maintenance services to ensure optimal performance and longevity of the device

Product training and education to ensure proper use and understanding of the device

Access to replacement parts and accessories to keep your nebulizer functioning at its best

24/7 customer support to assist with any questions or concerns you may have

We are committed to providing our customers with the highest level of support and service, so you can feel confident in your purchase of the Portable Nebulizer Machine.

Packing and Shipping:

Product Packaging:

1 Portable Nebulizer Machine

1 Mouthpiece

1 Adult Mask

1 Child Mask

1 USB Cable

1 User Manual

Shipping:

Shipping Method: Standard Shipping

Estimated Delivery Time: 3-5 business days

Shipping Cost: Free

Tracking Number: Provided

FAQ:

A: The brand name of this nebulizer machine is Janejoy.

Q: What is the model number of this nebulizer machine?

A: The model number of this nebulizer machine is J04.

Q: Where is this nebulizer machine made?

A: This nebulizer machine is made in China.

Q: Is this nebulizer machine portable?

A: Yes, this nebulizer machine is portable.

Q: What type of medication can be used with this nebulizer machine?

A: This nebulizer machine can be used with both liquid and powder medication.



Guangzhou Janejoy Medical Technology Co.,Ltd

+8613535187404

JaneJoy@therapy.org.cn

medicalcare-products.com

Building A2, Loteam Creative Park, No. 251, Kehua Street, Tianhe District, Guangzhou, China