



## Home Use Breathing Lung Expansion & Mucus Removal Device For Children And Adults

### Our Product Introduction

#### Basic Information

- Place of Origin: China
- Brand Name: JaneJoy
- Certification: FDA and CE
- Model Number: J-11
- Minimum Order Quantity: 72sets
- Price: 13USD
- Packaging Details: 72sets one carton packing
- Delivery Time: 25days lead time
- Payment Terms: T/T
- Supply Ability: 10000pcs per month



#### Product Specification

- Features: New Design
- Color: Blue And Green
- Usage: Breathing Exercises
- Resistance Levels: Adjustable
- Target Muscle Group: Respiratory Muscles
- Suitable For: All Ages
- Package Contents: 1 Breathing Muscle Trainer, Instruction Manual
- Material: ABS
- Size: One Size
- Warranty: 1 Year
- Benefits: Improved Breathing Efficiency
- Application: Rehab Center /Hospital /Home
- Function 1: Increases The Volume Of Lung Oxygen Intake.
- Function 2: Lesser Breath Holding Capacity



## Product Description

### Products Description

1-6 adjustable dials to practice inspiratory muscles and 1-5 levels of practicing expiratory muscles at one time. Good for Asthmatics, COPD, as well as breathing disordered patients.

<b>Product name</b>	Respiratory Muscle Trainer
<b>MOQ</b>	72 set
<b>Product structure</b>	Main body and blowpipe
<b>Use crowd</b>	Children, adults, the elderly
<b>Logo</b>	To be customized
<b>Color</b>	To be customized
<b>Sample</b>	Available
<b>Support OEM</b>	Support customized colors, packing, etc.

#### Respiratory Muscle Trainer

- Increases the volume of lung oxygen intake.
- Longer breath-holding capacity.
- Get rid of the trapped air in the alveoli for better gas exchange.
- Improve breathing and reduce shortness of breath.
- Help to improve and maintain lung hygiene by helping to remove pollutants from the lungs.
- Slowing down the build-up of lactate in the muscles.
- Protection box provided to keep the device from scratches, collision, dust, bacteria, and water.
- Low cost and effective option for mucus clearance.
- Short sessions required for maximum effectiveness.
- Keep the airways open.
- Healthy people of children and adults, respiratory disease sufferers, smokers, swimmers & deep sea divers, athletes, singers, and wind instrument musicians are suitable to use.

---

---

---

---

---

---

---

---

---

---



**Guangzhou Janejoy Medical Technology Co., Ltd**

+8613535187404

JaneJoy@therapy.org.cn

medicalcare-products.com

Building A2, Loteam Creative Park, No. 251, Kehua Street, Tianhe District, Guangzhou, China