



## Adjustable Resistance Lung Exerciser Device Wholesale Breathing Exercise Equipment Lung Expansion Training Sports Entert

### Our Product Introduction

for more products please visit us on [medicalcare-products.com](http://medicalcare-products.com)

#### Basic Information

- Place of Origin: China
- Brand Name: JaneJoy
- Certification: FDA and CE
- Model Number: J-15
- Minimum Order Quantity: 1000
- Price: 1.19USD one set
- Packaging Details: one set with color box, carton packing
- Delivery Time: 25days lead time
- Payment Terms: T/T
- Supply Ability: 10000pcs per month



#### Product Specification

- Size: One Size Fits All
- Color: Blue,white , Black Red
- Material: ABS
- Resistance Levels: Adjustable
- Suitable For: Athletes, Singers, Patients With Respiratory Conditions
- Weight: Lightweight
- Usage: Lung Breathing Exercises
- Warranty: 1 Year
- Benefits: Improves Lung Capacity, Strengthens Respiratory Muscles, Enhances Breathing Control
- Targeted Muscles: Diaphragm, Intercostal Muscles
- Includes: Instruction Manual, Carrying Case
- Highlight: **adjustable lung breathing trainer , adjustable breathing muscle trainer , adjustable resistance lung breathing trainer**



## Overview

### Products Description

Product name	Lung breathing trainer for respiratory muscle training
MOQ	1000 pieces
Product structure	Main body and blowpipe
Use crowd	Children,adults,the elderly
Logo	To be customized
Color	To be customized
Sample	Available
Support OEM	Support customized colors, packing, etc.

#### Principle of Breath Trainer

The breathing exerciser is used for resistance training through breathing holes of different sizes,which can enhance the strength and endurance of the respiratory muscles.Thereby reversing or reducing the state of respiratory muscle weakness and improving ventilation efficiency.Exhalation resistance training imitates lip-respiratory breathing,which can improve alveolar ventilation and gas exchange,help increase the body's alveolar ventilation and blood oxygen saturation,and improve lung function.

---

---

---



**Guangzhou Janejoy Medical Technology Co.,Ltd**



+8613535187404



JaneJoy@therapy.org.cn



medicalcare-products.com

Building A2, Loteam Creative Park, No. 251, Kehua Street, Tianhe District, Guangzhou, China