



Mini Ultrasonic Portable Mesh Nebuliser Inhaler Silent Asthma Atomizer For Kids

Our Product Introduction

For more products please visit us on medicalcare-products.com

Basic Information

- Place of Origin: China
- Brand Name: Janejoy
- Model Number: J04



Product Specification

- Type: Ultrasonic
- Power Supply Mode: USB Recharge
- Customer Logo: Customized
- Warranty: 1 Year
- Instrument Classification: Class II
- Material: ABS
- Atomization Rate: $\geq 0.2\text{ML/min}$
- Noise Level: Less Than 60dB
- Highlight: **Class II Ultrasonic Inhaler Mesh Nebulizer, 0.75w Ultrasonic Portable Mesh Nebuliser, USB Recharge Portable Inhaler Nebulizer**



More Images



Product Description

Portable Mini Hands Free Inhaler Silent Asthma Inhalator Atomizer Ultrasonic Mesh Nebulizer Medical Machine For Kids

Products Description

Name	Nebulizer
Model	U210
Dimensions	67*42*116mm
Power	0.75w
Power supply	USB charging
Atomization rate	>0.25ml/minute
Atomized median particle size	MMAD 4.0um
Weight	948g
Ambient temperature	10-40°C
Maximum water temperature in medicine cup	50°C Maximum
Work noise	50dB Maximum



Guangzhou Janejoy Medical Technology Co.,Ltd

+8613535187404

JaneJoy@therapy.org.cn

medicalcare-products.com

Building A2, Loteam Creative Park, No. 251, Kehua Street, Tianhe District, Guangzhou, China