



Electrical Stimulation Pelvic Floor Muscle Stimulator Vaginal Trainer With App Controlled

Our Product Introduction

for more products please visit us on medicalcare-products.com

Basic Information

- Place of Origin: China
- Brand Name: JaneJoy
- Certification: FDA and CE
- Model Number: J-20
- Minimum Order Quantity: 100
- Price: 33.4USD one set
- Packaging Details: one set one color box , carton packing
- Delivery Time: 25days lead time
- Payment Terms: , T/T
- Supply Ability: 10000pcs per month



Product Specification

- Material: ABS
- Warranty: 1year
- After-sale Service: Return And Replacement
- Instrument Classification: Class II
- APPLICATION: Pelvic Floor Trainer For Women Postpartum
- Packing: 1 Set/box
- Function: Muscle Trainer For Pelvic Floor Muscle Repair
- MOQ: 100sets
- Service: ODM OEM Service
- Color: White
- Highlight: **App Controlled Pelvic Floor Muscle Stimulator, Electrical Stimulation Pelvic Vaginal Trainer, Wireless Electrode Electronic Pelvic Floor Stimulator**



More Images



Product Description

Electrical Stimulation Biological Feedback Pelvic Floor Muscle Exercise Stimulator Vaginal Trainer With Optional App

Overview

Product Description

XFT-0010EB wireless pelvic muscle trainer adopts low-frequency electrical stimulation technology. Electrical stimulation of the intracavity electrodes that controlled by the Host and APP to instead of nerve impulses generated by the brain to produce equal tension or equal length contraction of the muscle strength training, thereby improving urinary incontinence and upgrading the quality of sex life.

Model No.:	XFT-0010EB
Shipment	By express/air/sea can be optional
Gift box size	19.5*16.5*7cm
Carton size	50*40*33 cm
Qty/CTN	30pcs
Weight/CTN	28kgs

The XFT-0010EB Wireless Pelvic Muscle Trainer mainly consists of a trainer, an intracavity electrode, a power cord and an APP.

Features •Minimalist dual probes design. •Scientifically equipped with an APP. •Portable self-contained carrying case. •Maintenance-free probe, can work for 30 hours. •Streamlined waterproof design of the probe, comfortable to use. •The APP supports to adjust the working parameters as required. •Wireless electrode, the trainer and the probe connect together via Bluetooth. •Using low-frequency electrical stimulation technology, which is a gentle electrical pulse.



Guangzhou Janejoy Medical Technology Co.,Ltd



+8613535187404



JaneJoy@therapy.org.cn



medicalcare-products.com

Building A2, Loteam Creative Park, No. 251, Kehua Street, Tianhe District, Guangzhou, China

