



Bladder Control Kegel Pelvic Trainer Pelvic Floor Stimulator For Incontinence

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: JaneJoy
- Certification: FDA and CE
- Model Number: J-23
- Minimum Order Quantity: 100
- Price: 90USD one set
- Packaging Details: one set one color box , carton packing
- Delivery Time: 25days lead time
- Payment Terms: , T/T
- Supply Ability: 10000pcs per month



Product Specification

- Material: ABS
- Warranty: 1year
- After-sale Service: Return And Replacement
- Instrument Classification: Class II
- APPLICATION: Pelvic Floor Trainer For Women Postpartum
- Packing: 1 Set/box
- Function: Muscle Trainer For Pelvic Floor Muscle Repair
- MOQ: 100sets
- Service: ODM OEM Service
- Color: White
- Highlight: **Bladder Control Kegel Pelvic Trainer,
Bladder Control Pelvic Floor Stimulator,
Custom Pelvic Floor Trainer For Incontinence**



for more products please visit us on medicalcare-products.com

Product Description

Best Seller Incontinence Stimulator Pelvic Floor Muscle Exercise Trainer /Pelvic Floor Stimulator /Pelvic Floor Exercise

Overview

Product Description

XFT-0010EA Pelvic Muscle Trainer mainly consists of a trainer, an intracavity electrode, a power cord and an APP software.

Model No.:	XFT-0010EA
Pulse width:	100-300 (±10%)
Frequency:	10-50Hz (±10%)
Output strength (load is 500Ω):	0-60mA
Training time:	6-30 min

No.	Parts	Quantity
1	Trainer	1 pcs
2	Intracavity electrode	1 pcs
3	APP software	1 pcs
4	Power Adapter	1pcs



Guangzhou Janejoy Medical Technology Co.,Ltd

+8613535187404

JaneJoy@therapy.org.cn

medicalcare-products.com

Building A2, Loteam Creative Park, No. 251, Kehua Street, Tianhe District, Guangzhou, China