



## Female Smart Wireless Pelvic Floor Muscle Trainer With App Controlled OEM ODM

Our Product Introduction

for more products please visit us on [medicalcare-products.com](http://medicalcare-products.com)

### Basic Information

- Place of Origin: China
- Brand Name: JaneJoy
- Certification: FDA and CE
- Model Number: J-23
- Minimum Order Quantity: 100
- Price: 90USD one set
- Packaging Details: one set one color box , carton packing
- Delivery Time: 25days lead time
- Payment Terms: , T/T
- Supply Ability: 10000pcs per month



### Product Specification

- Material: ABS
- Warranty: 1year
- After-sale Service: Return And Replacement
- Instrument Classification: Class II
- APPLICATION: Pelvic Floor Trainer For Women Postpartum
- Packing: 1 Set/box
- Function: Muscle Trainer For Pelvic Floor Muscle Repair
- MOQ: 100sets
- Service: ODM OEM Service
- Color: White
- Highlight: **App Connect Pelvic Floor Muscle Trainer,  
OEM Smart Pelvic Floor Trainer,  
ODM App Controlled Pelvic Floor Trainer**



## Product Description

### New Product Portable App Connect Wireless Female Improve Incontinence Pelvic Floor Muscle Tighten Trainer

#### Overview

#### Product Description

XFT-0010EA Pelvic Muscle Trainer mainly consists of a trainer, an intracavity electrode, a power cord and an APP software.

<b>Model No.:</b>	XFT-0010EA
<b>Pulse width:</b>	100-300 (±10%)
<b>Frequency:</b>	10-50Hz (±10%)
<b>Output strength (load is 500Ω):</b>	0-60mA
<b>Training time:</b>	6-30 min

No.	Parts	Quantity
1	Trainer	1 pcs
2	Intracavity electrode	1 pcs
3	APP software	1 pcs
4	Power Adapter	1 pcs



Guangzhou Janejoy Medical Technology Co.,Ltd

+8613535187404

JaneJoy@therapy.org.cn

medicalcare-products.com

Building A2, Loteam Creative Park, No. 251, Kehua Street, Tianhe District, Guangzhou, China