



## OEM Wired Physiotherapy Tens Electrodes Pad for Body Massager Muscle Stimulator

Our Product Introduction

For more products please visit us on [medicalcare-products.com](http://medicalcare-products.com)


### Basic Information

- Place of Origin: China
- Brand Name: Janejoy
- Certification: FDA and CE
- Model Number: J32
- Minimum Order Quantity: 5000
- Price: 0.17USD/pcs
- Packaging Details: 2000pcs one carton
- Delivery Time: 35days lead time



### Product Specification

- Latex-Free: Yes
- Reusable Time: Up To 30 Uses
- Flexible: Yes
- Easy To Apply: Yes
- Self-Adhesive: Yes
- Product: Physiotherapy Tens Electrode Pads Wire
- Gel Type: China ,Japan And USA Type Gel
- Adhesive: Gel
- Highlight: **OEM Tens Electrodes Pad ,  
Wired Physiotherapy Tens Electrodes Pad ,  
Soft Cloth Physiotherapy Electrode Pads**



### More Images



## Product Description

### Relax Body Massager Electrical Muscle Stimulator Wired Physiotherapy Tens Electrodes Pads

#### Product Features:

1. Soft cloth, foam, tan tricot or other colours backing for selection
2. Durable structure allows patients removing by pulling up from the wire without damaging the pad
3. Pigtail and snap connection (optional)
4. Elastic socket fit male pin size from 1.95 to 2.15mm
5. Water base gel makes it skin friendly and allergy free
6. Soft cloth backing conform to any area of body
7. Pads, liners and pigtails can be printed your logo
8. Various sizes and any shapes are available
9. Self-adhesive, reusable, bio-compatible, No hot spot & over saturation
10. Various connector: round molded/heat shrink/eziehold flat/moulded with transparent lead wire
11. CE certified LATEX-FREE
12. Shelf lifespan: at least 2 years

#### Products Description

<b>Sample</b>	Acceptable
<b>Function</b>	Reduce Pain,Physical therapy
<b>Part Used</b>	Body , Arm , leg , foot , neck , back , waist etc.
<b>Application</b>	For TENS unit, for EMS unit,for medical device
<b>OEM ODM service</b>	Support custom LOGO , Packaging, color, shape
<b>Delivery time</b>	About 45 days after order confirmed
<b>Payment terms</b>	T/T,L/C, Western Union, Paypal, Bank account , Alibaba Trade Assurance ( Visa, Mastercard, e-Checking, T/T)



Guangzhou Janejoy Medical Technology Co.,Ltd



+8613535187404



JaneJoy@therapy.org.cn



medicalcare-products.com

Building A2, Loteam Creative Park, No. 251, Kehua Street, Tianhe District, Guangzhou, China