




Reusable Tens Electrode Pad For Relieve Pain Body Massage Physical Therapy OEM ODM

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: Janejoy
- Certification: FDA and CE
- Model Number: J35
- Minimum Order Quantity: 5000
- Price: 0.2USD/pcs
- Packaging Details: 1000pcs one carton
- Delivery Time: 35days lead time



Product Specification

- Latex-Free: Yes
- Reusable Time: Up To 30 Uses
- Flexible: Yes
- Easy To Apply: Yes
- Self-Adhesive: Yes
- Product: Physiotherapy Tens Electrode Pads Wire
- Gel Type: China ,Japan And USA Type Gel
- Adhesive: Gel
- Material: Soft Cloth
- Color: Brown
- Highlight: **OEM Tens Electrode Pad ,
Body Massage Tens Electrode Pads,
Physical Therapy Tens Electrode Pads**



Product Description

Economical Body Massage Physical Therapy Reusable Relieve Pain Tens Electrode Pads

Products Description

Sample	Acceptable
Function	Reduce Pain,Physical therapy
Part Used	Body , Arm , leg , foot , neck , back , waist etc.
Application	For TENS unit, for EMS unit,for medical device
OEM ODM service	Support custom LOGO , Packaging, color, shape
Delivery time	About 45 days after order confirmed
Payment terms	T/T,L/C, Western Union, Paypal, Bank account , Alibaba Trade Assurance (Visa, Mastercard, e-Checking, T/T)

Product Features: 1. Soft cloth, foam, tan tricot or other colours backing for selection 2. Durable structure allows patients removing by pulling up from the wire without damaging the pad 3. Pigtail and snap connection (optional) 4. Elastic socket fit male pin size from 1.95 to 2.15mm 5. Water base gel makes it skin friendly and allergy free 6. Soft cloth backing conform to any area of body 7. Pads, liners and pigtails can be printed your logo 8. Various sizes and any shapes are available 9. Self-adhesive, reusable, bio-compatible, No hot spot & over saturation 10. Various connector: round molded/heat shrink/eziehold flat/moulded with transparent lead wire 11. CE certified LATEX-FREE 12. Shelf lifespan: at least 2 years



Guangzhou Janejoy Medical Technology Co.,Ltd



+8613535187404



JaneJoy@therapy.org.cn



medicalcare-products.com

Building A2, Loteam Creative Park, No. 251, Kehua Street, Tianhe District, Guangzhou, China