




## Class I Physical Therapy Equipments Adhesive Electrode Pads For EMS / TENS Massager

### Our Product Introduction

#### Basic Information

- Place of Origin: China
- Brand Name: Janejoy
- Certification: FDA and CE
- Model Number: J35
- Minimum Order Quantity: 5000
- Price: 0.2USD/pcs
- Packaging Details: 1000pcs one carton
- Delivery Time: 35days lead time



#### Product Specification

- Latex-Free: Yes
- Reusable Time: Up To 30 Uses
- Flexible: Yes
- Easy To Apply: Yes
- Self-Adhesive: Yes
- Product: Physiotherapy Tens Electrode Pads Wire
- Gel Type: China ,Japan And USA Type Gel
- Adhesive: Gel
- Material: Soft Cloth
- Color: Brown
- Highlight: **TENS Massager Adhesive Electrode Pads,  
EMS Massager Adhesive Electrode Pads,  
Physical Therapy Equipments Electrode Pads**

## Product Description

### Custom Class I Physical Therapy Equipments Adhesive Electrode Pad For EMS/ TENS Massager

#### Products Description

<b>Sample</b>	Acceptable
<b>Function</b>	Reduce Pain,Physical therapy
<b>Part Used</b>	Body , Arm , leg , foot , neck , back , waist etc.
<b>Application</b>	For TENS unit, for EMS unit,for medical device
<b>OEM ODM service</b>	Support custom LOGO , Packaging, color, shape
<b>Delivery time</b>	About 45 days after order confirmed
<b>Payment terms</b>	T/T,L/C, Western Union, Paypal, Bank account , Alibaba Trade Assurance ( Visa, Mastercard, e-Checking, T/T)

#### Product Features:

1. Soft cloth, foam, tan tricot or other colours backing for selection
2. Durable structure allows patients removing by pulling up from the wire without damaging the pad
3. Pigtail and snap connection (optional)
4. Elastic socket fit male pin size from 1.95 to 2.15mm
5. Water base gel makes it skin friendly and allergy free
6. Soft cloth backing conform to any area of body
7. Pads, liners and pigtails can be printed your logo
8. Various sizes and any shapes are available
9. Self-adhesive, reusable, bio-compatible, No hot spot & over saturation
10. Various connector: round molded/heat shrink/eziehold flat/moulded with transparent lead wire
11. CE certified LATEX-FREE
12. Shelf lifespan: at least 2 years



**Guangzhou Janejoy Medical Technology Co.,Ltd**

+8613535187404

JaneJoy@therapy.org.cn

medicalcare-products.com

Building A2, Loteam Creative Park, No. 251, Kehua Street, Tianhe District, Guangzhou, China