



5cm*5cm Reusable Tens Electrode Pad , Soft Cloth Tens Machine Electrode Pads

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: JaneJoy
- Certification: FDA and CE
- Model Number: J-36
- Minimum Order Quantity: 5000
- Price: 0.17USD one pcs
- Packaging Details: 4pcs one pack , 2000pcs one carton
- Delivery Time: 40days lead time
- Supply Ability: 10000000pcs per month



Product Specification

- Material: Soft Cloth
- Color: Custom
- Size: 5*5cm
- Type: Pre-Wired Type
- OEM Service: Support Custom LOGO , Packaging, Color, Shape
- Packing: 4pcs/pack In A Resealable Zip Lock Bag
- Hydrogel: Biocompatible Conductive Hydrogel
- Shape: Square/custom
- Logo Printed: The Paper Insert, Plastic Bag, Electrode's Backing
- Application: For TENS Unit, For EMS Unit,for Medical Device
- Function: Reduce Pain
- Connector: Round Moulded Connector (other Backing, Conductor And Connector Optional)



Product Description

5cm*5cm Reusable Self-Adhering Tens Unit Electrode Pads For Ems Electric Muscle Stimulator

Key attributes

Industry-specific attributes

Warranty	3 years
After-sale Service	Return and Replacement

Other attributes

Place of Origin	Guangdong, China
Instrument classification	Class I
Brand Name	Bytech
Model Number	CWN2505(custom logo)
TENS electrode pad	50*50mm
Backing	White fabric
Product name	tens replacement pad
Conductive layer	Carbon film
Connector	Pre-wire
Application	for TENS/EMS/NMES/IF/therapy devices
Quality approval	CE/ISO
Warranty	3 years
MOQ	5000pcs
Service	ODM OEM Service

Packaging and delivery

Packaging Details	4pcs/pack into a zip bag, but customer's packing available
Port	Guangzhou

Supply Ability

Supply Ability	2000000 Piece/Pieces per Month tens replacement pad
----------------	---

[Show less](#)

Lead time

Quantity (pieces)	1 - 1000	> 1000
Lead time (days)	15	To be negotiated

Customization

Customized logo

Min. order: 10000

Customized...
packaging

Min. order: 10000

Graphic...
customization


Min. order: 10000

For more customization details, [message supplier](#)

Overview



Product Description



Product name	Electrodes Pads for TENS/EMS
Model	CWN2505
Size	50*50mm; 2"x2"
Shape	Square/custom
Backing	White fabric/Soft cloth/custom
Conductor	Super conductive carbon film
Connector	Round moulded connector (other backing, conductor and connector optional)
Type	Pre-Wired type
Hydrogel	Biocompatible conductive hydrogel
Logo printed	The paper insert, plastic bag, electrode's backing
Function	Reduce Pain
Part Used	Body , Arm , leg , foot , neck , back , waist etc.
Application	For TENS unit, for EMS unit,for medical device
Packing	4pcs/pack in a resealable zip lock bag
OEM ODM service	Support custom LOGO , Packaging, color, shape

Product Features

1. Special for Ac only 2. Reusable and self-adhering 3. Safe and comfortable 4. Durable 5. Stable 6. Variety of shapes and sizes



Product instructions:

1. Please clean the surface of the skin before use to remove surface dirt, sweat and dust. 2. Correctly connect the massage patch to the therapy instrument wire, remove the protective film of the massage patch, paste it on the massage site, and start the routine process. 3. After each massage, please turn off the machine first, then remove the wire plug, remove the massage patch, paste the protective film, and store it properly for next use. 4. The general service life of massage stickers is 100-200 times. (The service life of the product is related to skin characteristics and usage methods) 5. Please seal and store

in a cool place

Recommended by seller

[view more](#)

>>

New style practical tens unit self adhesive reusable physical therapy electrode pad
€0.184 - €0.9199 / piece
5000 pieces

Body stimulation physical therapy replacement massage self adhesive tens unit electrode pads
€0.092 - €0.9199 / piece
5000 pieces

CE certified rehabilicare product re-usable self adhesive tens unit massage electrode pads
€0.092 - €0.9199 / piece
5000 pieces

Custom round non-woven TENS self-adhesive physiotherapy snap electrodes patches pad
€0.46 - €0.9199 / pack
5000 packs

We have a professional R&D and design team
Welcome OEM ODM Service
Support Custom color, LOGO, packaging, material...

Application



Guangzhou Janejoy Medical Technology Co.,Ltd



+8613535187404



JaneJoy@therapy.org.cn



medicalcare-products.com

Building A2, Loteam Creative Park, No. 251, Kehua Street, Tianhe District, Guangzhou, China