



## 6x8cm Tens Unit Replacement Reusable Electrode Pads For Electronic Pulse Massager

### Our Product Introduction

#### Basic Information

- Place of Origin: China
- Brand Name: JaneJoy
- Certification: FDA and CE
- Model Number: J-36
- Minimum Order Quantity: 5000
- Price: 0.46USD one pcs
- Packaging Details: 4pcs one pack , 2000pcs one carton
- Delivery Time: 40days lead time
- Supply Ability: 10000000pcs per month



#### Product Specification

- Material: Soft Cloth
- Color: Custom
- Size: 6\*8cm
- Type: Pre-Wired Type
- OEM Service: Support Custom LOGO , Packaging, Color, Shape
- Packing: 4pcs/pack In A Resealable Zip Lock Bag
- Hydrogel: Biocompatible Conductive Hydrogel
- Shape: Square/custom
- Logo Printed: The Paper Insert, Plastic Bag, Electrode's Backing
- Application: For TENS Unit, For EMS Unit,for Medical Device
- Function: Reduce Pain
- Connector: Round Moulded Connector (other Backing, Conductor And Connector Optional)

## Product Description

### 6\*8cm Tens Unit Replacement Reusable Electrode Pads For Electronic Pulse Massager

#### Products Description

<b>Product name</b>	Electrode pads
<b>Size</b>	60*80mm or coustom
<b>Sample</b>	Acceptable
<b>Function</b>	Reduce Pain,Physical therapy
<b>Part Used</b>	Body , Arm , leg , foot , neck , back , waist etc.
<b>Application</b>	For TENS unit, for EMS unit,for medical device
<b>OEM ODM service</b>	Support custom LOGO , Packaging, color, shape
<b>Delivery time</b>	About 45 days after order confirmed
<b>Payment terms</b>	T/T,L/C, Western Union, Paypal, Bank account , Alibaba Trade Assurance ( Visa, Mastercard, e-Checking, T/T)



Guangzhou Janejoy Medical Technology Co.,Ltd



+8613535187404



JaneJoy@therapy.org.cn



medicalcare-products.com

Building A2, Loteam Creative Park, No. 251, Kehua Street, Tianhe District, Guangzhou, China