



60x80mm Physiotherapy Electrode Pads For Tens Unit , Tens EMS Pads Custom Logo

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: JaneJoy
- Certification: FDA and CE
- Model Number: J-36
- Minimum Order Quantity: 5000
- Price: 0.46USD one pcs
- Packaging Details: 4pcs one pack , 2000pcs one carton
- Delivery Time: 40days lead time
- Supply Ability: 10000000pcs per month



Product Specification

- Material: Soft Cloth
- Color: Custom
- Size: 6*8cm
- Type: Pre-Wired Type
- OEM Service: Support Custom LOGO , Packaging, Color, Shape
- Packing: 4pcs/pack In A Resealable Zip Lock Bag
- Hydrogel: Biocompatible Conductive Hydrogel
- Shape: Square/custom
- Logo Printed: The Paper Insert, Plastic Bag, Electrode's Backing
- Application: For TENS Unit, For EMS Unit,for Medical Device
- Function: Reduce Pain
- Connector: Round Moulded Connector (other Backing, Conductor And Connector Optional)

Product Description

Premium Physiotherapy Self Adhesive Tens Ems Electrical Electrode Pads

Products Description

Product name	Electrode pads
Size	60*80mm or coustom
Sample	Acceptable
Function	Reduce Pain,Physical therapy
Part Used	Body , Arm , leg , foot , neck , back , waist etc.
Application	For TENS unit, for EMS unit,for medical device
OEM ODM service	Support custom LOGO , Packaging, color, shape
Delivery time	About 45 days after order confirmed
Payment terms	T/T,L/C, Western Union, Paypal, Bank account , Alibaba Trade Assurance (Visa, Mastercard, e-Checking, T/T)



Guangzhou Janejoy Medical Technology Co.,Ltd



+8613535187404



JaneJoy@therapy.org.cn



medicalcare-products.com

Building A2, Loteam Creative Park, No. 251, Kehua Street, Tianhe District, Guangzhou, China