




## Facial Massager Non Woven Adhesive Electrode Pad For Face Beauty Device

Our Product Introduction

for more products please visit us on [medicalcare-products.com](http://medicalcare-products.com)


### Basic Information

- Place of Origin: China
- Brand Name: JaneJoy
- Certification: FDA and CE
- Model Number: J-36
- Minimum Order Quantity: 5000
- Price: 1USD one pcs
- Packaging Details: 4pcs one pack , 2000pcs one carton
- Delivery Time: 40days lead time
- Supply Ability: 10000000pcs per month



### Product Specification

- Material: Soft Cloth
- Color: Custom
- Type: Pre-Wired Type
- OEM Service: Support Custom LOGO , Packaging, Color, Shape
- Packing: 4pcs/pack In A Resealable Zip Lock Bag
- Hydrogel: Biocompatible Conductive Hydrogel
- Shape: Square/custom
- Logo Printed: The Paper Insert, Plastic Bag, Electrode's Backing
- Application: For TENS Unit, For EMS Unit, for Medical Device
- Function: Reduce Pain
- Connector: Round Moulded Connector (other Backing, Conductor And Connector Optional)
- Backing: White Fabric/Soft Cloth/custom



### More Images



## Product Description

### Facial Massager Non Woven Adhesive Electrode Pad For Face Beauty Device

#### Products Description

<b>Product name</b>	face electrode pads
<b>Sample</b>	Acceptable
<b>Function</b>	Reduce Pain,Physical therapy
<b>Part Used</b>	Face
<b>Application</b>	For TENS unit, for EMS unit,for medical device
<b>OEM ODM service</b>	Support custom LOGO , Packaging, color, shape
<b>Delivery time</b>	About 45 days after order confirmed
<b>Payment terms</b>	T/T,L/C, Western Union, Paypal, Bank account , Alibaba Trade Assurance ( Visa, Mastercard, e-Checking, T/T)

**Product Features:** 1. Soft cloth, foam, tan tricot or other colours backing for selection 2. Durable structure allows patients removing by pulling up from the wire without damaging the pad 3. Pigtail and snap connection (optional) 4. Elastic socket fit male pin size from 1.95 to 2.15mm 5. Water base gel makes it skin friendly and allergy free 6. Soft cloth backing conform to any area of body 7. Pads, liners and pigtails can be printed your logo 8. Various sizes and any shapes are available 9. Self-adhesive, reusable, bio-compatible, No hot spot & over saturation 10. Various connector: round molded/heat shrink/eziehold flat/moulded with transparent lead wire 11. CE certified LATEX-FREE 12. Shelf lifespan: at least 2 years



**Guangzhou Janejoy Medical Technology Co.,Ltd**



+8613535187404



JaneJoy@therapy.org.cn



medicalcare-products.com

Building A2, Loteam Creative Park, No. 251, Kehua Street, Tianhe District, Guangzhou, China