



CE Portable Mesh Nebulizer With Nebulization Adjustable Inhaler Medication Nebulizer Machine Water Temperature $\leq 50^{\circ}\text{C}$

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: JaneJoy
- Certification: CE and FDA
- Model Number: UN200
- Minimum Order Quantity: 200
- Price: 8
- Packaging Details: 50sets one carton
- Delivery Time: 25days lead time
- Payment Terms: T/T
- Supply Ability: 1000000pcs per month



Product Specification

- Material: ABS
- Color: White
- Power: $<2\text{W}$
- Nebulization Rate: $\geq 0.25\text{ ML/min}$
- Respirable Fraction: $>60\%$
- Vibration Frequency: $110\text{kHz} \pm 10\%$
- Volume: Max.10mL
- Water Temperature: ≤ 50
- Warranty: 1 Year
- Size: $67 \times 45 \times 105\text{mm}$
- Product Weight: $95\text{g} \pm 10\%$ (Without Batteries)
- Suitable For: Adults And Children
- Noise Level: Less Than 50 DB
- Voltage: DC3V USB 5V OR DC3.7V Rechargeable Lithium Battery

Product Description

Our Product Introduction

for more products please visit us on medicalcare-products.com

Products Description

Model	UN200
Voltage	DC3V USB 5V OR DC3.7V Rechargeable Lithium Battery
Power	<2W
Nebulization rate	≥0.25 mL/min
Respirable fraction	>60%
Vibration frequency	110kHz±10%
Volume	Max.10mL
Water Temperature	≤50
Product weight	95g±10%(Without batteries)
Size (L×W×H)	67×45×105mm

Shipping & payment



Guangzhou Janejoy Medical Technology Co.,Ltd



+8613535187404



JaneJoy@therapy.org.cn



medicalcare-products.com

Building A2, Loteam Creative Park, No. 251, Kehua Street, Tianhe District, Guangzhou, China