



Household Automatic Drug Compressor Mesh Nebulizer with High Nebulization Rate of ≥ 0.25 mL/min

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: JaneJoy
- Certification: CE and FDA
- Model Number: UN100
- Minimum Order Quantity: 200
- Price: 8
- Packaging Details: 50sets one carton
- Delivery Time: 25days lead time
- Payment Terms: T/T
- Supply Ability: 1000000pcs per month



Product Specification

- Material: ABS
- Color: White
- Power: <2W
- Nebulization Rate: ≥ 0.25 ML/min
- Respirable Fraction: >60%
- Vibration Frequency: 110kHz \pm 10%
- Volume: Max.10mL
- Water Temperature: ≤ 50
- Warranty: 1 Year
- Size: 67 \times 45 \times 105mm
- Product Weight: 95g \pm 10%(Without Batteries)
- Suitable For: Adults And Children
- Noise Level: Less Than 50 DB
- Voltage: DC3V USB 5V OR DC3.7V Rechargeable Lithium Battery



for more products please visit us on medicalcare-products.com


Product Description




Our Product Introduction

for more products please visit us on medicalcare-products.com

Products Description

Model	UN100
Voltage	DC 3V USB 5V OR DC3.7V Rechargeable Lithium Battery
Power	<2W
Particle size	MMAD:
Nebulization rate	≥0.25 mL/min
Particle distribution(%)	(0.5-5μm)
Respirable fraction	>60%
Vibration frequency	110kHz±10%
Volume	Max.10mL
Water Temperature	≤50
Working voice	≤50dB(A)


Guangzhou Janejoy Medical Technology Co.,Ltd

 +8613535187404
  JaneJoy@therapy.org.cn
  medicalcare-products.com
 Building A2, Loteam Creative Park, No. 251, Kehua Street, Tianhe District, Guangzhou, China